



# Intraoral Video Camera System

*Part Number D8000*

## Operator's Manual



**AIR  
TECHNIQUES**  
INC.



***FOREWORD***

---

Air Techniques has prepared this document as a guide to the proper use of the Astra Intraoral Video Camera System. Review and follow the guidelines included in this Operator’s Manual to ensure that your Astra System gives you the highest level of service. If any difficulties are encountered with this product, please contact your authorized dealer.

***TABLE OF CONTENTS***

---

	<u>Page</u>
General Safety . . . . .	3
Congratulations . . . . .	4
Important Information . . . . .	5
Key Parts Identification . . . . .	6
Component Description . . . . .	7
Unpacking and Inspection . . . . .	7
Pre-Installation Information . . . . .	8
Handpiece Functions . . . . .	9
System Setup . . . . .	10
Camera Portability . . . . .	12
Pre-Operation . . . . .	13
System Operation . . . . .	14
System Check . . . . .	15
Maintenance . . . . .	16
Specifications . . . . .	17
Physical Characteristics . . . . .	17
Accessories . . . . .	18
VISIX Imaging Software . . . . .	18
Warranty . . . . .	19
On-Line Warranty Registration . . . . .	19

This equipment has been designed to minimize exposure of personnel to hazards. While the Astra System is designed for safe operation, certain precautions must be observed. Use of the Astra System not in conformance with the instructions specified in this manual may result in permanent failure of the unit.

### **KNOWLEDGE OF WARNINGS AND CAUTIONS**

Users must exercise every precaution to ensure personnel safety, and be familiar with the warnings and cautions presented throughout this manual and summarized below.

#### **Warnings.**

- ☐ To prevent fire or electrical shock, do not expose this equipment to rain or moisture.
- ☐ Do not spray cleaning liquids or disinfectants directly on the camera handpiece or the Handpiece Cable. Use care not to allow liquids to run into internal circuitry.
- ☐ Do not wipe the surfaces using benzine, thinner, etc. as this may degrade the finish.

#### **Cautions.**

- ☐ There are no user-serviceable parts inside. Servicing should be performed by qualified dealer service representatives only.
- ☐ Connecting any device to the Astra System that does not meet the equivalent safety requirements of the system may reduce the safety effectiveness of the Astra System.
- ☐ Unplug the system components from power before performing cleaning procedures.

#### **Contraindications.**

- ☐ None known.

#### **Use of Accessory Equipment.**

All peripheral equipment (computers, monitors, printers, etc) connected to this equipment must comply with IEC 60601-1 and the resulting system must comply with IEC 60601-1 and/or IEC 60601-1-1 harmonized national standards when set-up within the immediate patient vicinity (ie. within 6 feet).

Accessory equipment outside the patient vicinity must comply with IEC 60601-1 or IEC 60950-1 as allowed by IEC 60601-1-1.

Use of ACCESSORIES or cables other than those specified or provided by Air Techniques may result in increased EMISSIONS or decreased IMMUNITY of the equipment.

#### **Do Not Attempt Internal Service.**

The interior of the camera is only accessible by removing hardware with tools and should only be opened and serviced by an authorized service technician.

Contact your local Air Techniques authorized dealer for service. Failure to heed this directive may result in equipment damage and voiding the warranty.

## GENERAL SAFETY

---

### Markings.

The following terms or symbols are used on the equipment or in this manual to denote information of special importance:



Alerts users to important Operating and Maintenance instructions. Read included documents.



Indicates type BF equipment in accordance with UL/IEC 60601-1



Alerts users that the camera cover must be used for only one patient and disposed of properly in accordance with local code.



**Medical Electrical Equipment**  
with respect to electrical shock, fire,  
mechanical and other specified hazards  
only in accordance with  
UL-60601-1, CAN/CSA C22.2 No.601.1  
66CA



Identifies the  
name of the  
manufacturer.



Indicates the  
month and year  
of manufacture



Indicates that the Astra complies with the Medical  
Device Directive 93/42/EEC.



Identifies the name of  
the manufacturer.



Authorized Representative for Medical  
Device Directive



Indicates date of manu-  
facture

## CONGRATULATIONS

---

Congratulations on your purchase of the Astra Intraoral Video Camera System, the latest camera in the dental video imaging product line from Air Techniques, a leading manufacturer of dental equipment since 1962.

The Astra Camera System design offers the flexible benefits of USB Plug-N-Play technology that is easily integrated into any practice. The very lightweight one-piece Astra Camera Handpiece is an extremely maneuverable instrument, which produces crisp, clear images. The Astra Intraoral Video Camera System is hereafter referred to as the Astra System in this manual.

Review and follow the guidelines included in this Operator's Manual to ensure that your Astra System gives you the highest level of service.

For product support and information on how to expand your Astra System, contact your authorized Air Techniques dealer.

### **General Notes.**

- ☐ All instructions in this manual form an integral part of the unit. They must be kept close to the unit and in readiness whenever required. Precise observance of these instructions is a pre-condition for use of the unit for the intended purpose and for its correct operation. This manual should be passed on to any future purchaser or operator.
- ☐ Safety of the operator as well as trouble-free operation of the unit are only ensured if use is made of original equipment parts. Moreover, use may only be made of those accessories that are specified in the technical documentation or that have been expressly approved and released by the manufacturer for the intended purpose. The manufacturer cannot warranty for the safety or proper functioning of this unit in the case where parts or accessories are used that are not supplied by the manufacturer.
- ☐ There is no guarantee against damage arising where parts or accessories are used that are not supplied by the manufacturer.
- ☐ Observe the usage and storage conditions.
- ☐ The Astra System may only be operated using authorized software.
- ☐ The manufacturer regard themselves as being responsible for the equipment with regard to safety, reliability and proper functioning only if assembly, resetting, changes or modifications and repairs have been carried out by an authorized dealer and if the equipment is used in conformity with the instructions contained in this manual.
- ☐ The device conforms to the relevant safety standards valid at this time.
- ☐ Any reprinting of the technical documentation, in whole or in part, is subject to prior written approval by the manufacturer.







### **Incorrect Usage**

- ☐ Any use that is not described in this manual as correct usage is considered as incorrect usage. The manufacturer is not to be held liable for any damage caused as a result of incorrect usage. The operator bears all risks.

**KEY PARTS IDENTIFICATION**

The Astra System consists of a Camera Handpiece, a Camera Handpiece Cable and associated accessories. The handpiece, cable, accessories and software may be purchased together under part number D8000. Some of the items may be purchased individually by the part numbers listed on page 18.

Astra System, Part Number D8000, includes:

Item	Qty	
Camera Handpiece	1	
Disposable Camera Covers Sample Pack for Camera Handpiece	50	
Handpiece Cable, 10 feet	1	
Handpiece Holder with Chair-side Adapter	1	
CD containing - PDF version of Operator's Manual USB Device Drivers Demonstration Program and Utilities Adobe Acrobat Reader	1	
Quick Start Instructions	1	

---

## ***COMPONENT DESCRIPTION***

The components of the Astra System are described below.

**Camera Handpiece** - The Camera Handpiece is extremely lightweight with a high resolution, high sensitivity, auto-exposure controlled CCD sensor and a high performance lens system illuminated by ultra-bright white LED lamps. The fixed focus lens has a broad depth of field enabling the camera to finely detail a section of a single tooth to a tour of the mouth. The Camera Handpiece has 2 keypad buttons that provide fingertip operating control of the system. The specific functions performed by the keypad buttons depend on the user-supplied compatible streaming video capture software application used.

**Camera Covers Sample Pack** - Provides a sample quantity of 20 disposable camera covers used as an effective barrier preventing any infection hazard to the patient.

**Handpiece Cable** - A 10 foot cable that is used to connect the Astra handpiece to the user's computer. One end of the cable has a keyed, 5-pin quick disconnect connector for connection to the Astra handpiece and the other end has a Type A USB connector for connection to a high speed, high power (500 mA) USB 2.0 port on the user's computer. This cable carries all power, video and keypad button control signals for the handpiece.

**Handpiece Holder** - The Handpiece Holder can be fastened in a convenient location for safe storage of the handpiece when it is not in use. In addition to storage, the handpiece holder also has built-in magnets that turn off the Astra Handpiece power. This feature conserves energy and equipment wear, allowing the system to be ready and operational when needed.

**Compact Disk** - A compact disk (CD) included with the Astra System that includes:

1. USB Device Drivers, Demonstration Program and Utilities
2. PDF version of Operator's Manual
3. Adobe Acrobat Reader

---

## ***UNPACKING AND INSPECTION***

Unpack each component of the Astra System and inspect for physical damage such as scratches, damaged connectors, etc. If any damage is noted, immediately notify your Air Techniques authorized dealer immediately so corrective action can be taken.

Make sure to save all packaging material in case repackaging and shipment is necessary.

Verify that all listed items were received. If any item is missing, notify your Air Techniques authorized dealer.

# PRE-INSTALLATION INFORMATION

**IMPORTANT:** When operating the Astra System connected to the high power (500 mA) USB 2.0 port on a Computer System, the computer must also be loaded with Air Techniques Authorized Twain or DirectX 9 compliant streaming video software application such as VISIX. Contact your dealer for available Computer Systems and software options.

**NOTE:** VISIX Imaging Software is fully compatible with Astra Systems.

**Minimum Computer System Requirements** - The Computer System (laptop or computer, monitor, etc.) and any related peripheral or other equipment, supplied by the user, or a third party, must comply with the requirements for accessory Equipment as specified in the appropriate IEC 60601-1 and/or IEC 60601-1-1 harmonized national standard. Unless otherwise stated, the components of the computer system must comply with the minimum requirements listed below.

- CPU Speed: Pentium-4, 2 GHz or higher
- Operating System: Windows XP Professional with Service Pack 3 or later for an Intel 32-bit processor;  
Microsoft Windows XP Professional 64-bit Edition with Service Pack 2 or later for an Intel 64-bit extended (x64) processor;  
Microsoft Windows Vista Business, Enterprise, or Ultimate with Service Pack 1 or later for an Intel 32-bit or an Intel 64-bit extended (x64) processor;  
Microsoft Windows 7 Professional, Enterprise, or Ultimate for an Intel 32-bit or an Intel 64-bit extended (x64) processor.
- System RAM: 256 MB
- Optical Drive: CD-ROM Capable
- Monitor : 800 x 600 resolution or higher
- Video Display Adapter: 16 MB video card with 800 x 600 pixel resolution and 32-bit color
- USB Port: USB 2.0 High-Speed, High-Power Port

## System Properties.

If unsure of the operating system version installed, check that it meets the necessary requirements by checking the **System Properties** window.

This is done simply by right clicking the My Computer icon. Selecting Properties from the menu list displays the **System Properties** window as shown. The installed operating system version is listed under the **General Tab**.

The **System Properties** window can also be displayed from the Desktop Start button. Just press the Start button and select **Settings→Control Panel** and then System.





### Astra Camera Handpiece



**Keypad Buttons** - Each keypad button provides the Astra Camera Handpiece fingertip control over the PC streaming video capture software via the Handpiece Cable. The specific functions performed by the keypad buttons depend on a software installation that supports the 2-button handpiece. Typical keypad actions provided below.

**NOTE:** Simultaneously depressing both keypad buttons for approximately 3 seconds deactivates the Camera Handpiece LED light source. This is commonly used when imaging an X-ray on a light box.

- |                          |  |
|--------------------------|--|
| (1) Top Keypad Button    | Press and release to send a <i>Button 1</i> command to the PC. While the specific function performed depends on the setting of the video capture software application installed on the PC, this command typically freezes and unfreezes the image that is displayed on the computer monitor. |
| (2) Bottom Keypad Button | Press and release to send a <i>Button 2</i> command to the PC. While the specific action depends on the video capture software installed, it typically saves the image displayed on the computer monitor to the computer hard disk.  |
| (3) Handpiece Connector  | A 5-pin connector socket that accepts connection of the keyed 5-pin quick disconnect plug end of the Handpiece Cable.  |

Installing the Astra System is as simple as making the necessary connections of the Handpiece Cable between a PC loaded with an Air Techniques Authorized DirectX 9 compliant user-supplied streaming video software application such as VISIX and the Camera Handpiece. Perform the following procedures to set up the Astra System.

**Notes:** Authorized Imaging Software, such as VISIX, supplied by the dealer or other company, must be installed on the computer in order to operate the Astra System.


A Microsoft® **Driver Not Signed** statement may appear. If it does, click **Continue** to proceed with installation.

**USB Device Driver Installation** - Before connecting the Astra System to your computer or attempting to use it for the first time, run the Setup program on the Astra System Drivers and Utilities Disk included with the Astra System. The supplied CD provides the device drivers needed to communicate between the Astra System and the associated user-supplied USB streaming video software application installed on the PC computer system. This must be completed before the Astra System will operate properly. Normally, this program runs automatically when the CD is inserted into the drive for the first time. If not, run the setup program located in the root directory of the CD (typically **D:\Astra.exe**).

Install the Standard Device Drivers by selecting Standard (default option) from the Setup program menu. If full device functionality is not present after the Standard Device Drivers are installed, the Legacy Device Drivers may need to be installed. More information can be found in the Installation Instructions and Notes file on the Astra System Drivers and Utilities Disk included with the system.

**Camera Handpiece Connection** - Carefully connect the Astra Camera Handpiece to the supplied Handpiece Cable as follows:


1. Using the gray molded connector end of the Handpiece Cable, align the connector key with the keyway of the Handpiece connector.
2. Insert straight into Handpiece connector until it securely snaps into place.



Camera Handpiece Connection

**USB Interface Connection** - Connect the USB Type A connector side of the Handpiece Cable to any available USB 2.0 High Speed, high power port on the user supplied Computer System. Always make sure to use the same USB port whenever re-connection of the USB cable is necessary.


**Important:** Use a High Power (500 mA) USB 2.0 Port Only



### **USB Interface Connection**

**Handpiece Holder with Chair-side Adapter** - Installing the Chair-side Adapter is as easy as deciding on a location for maximum convenience and securing it to an empty customer-supplied standard Handpiece Holder. Refer to the illustration below and install the Chair-side Adapter as follows:

1. Insert Handpiece Holder Adapter (1) into the wide end of Grip (2).
2. Press Grip (2) into the top of the chosen Handpiece Holder.
3. Insert Expander (3) into the bottom of Grip (2)
4. Insert and push the Joint Connector Hex Bolt (4) into Expander (3) engaging the thread of Handpiece Holder Adapter (1) by turning the Joint Connector Hex Bolt (4) clockwise.
5. Using the supplied Short Arm Allen Wrench (5), tighten Joint Connector Hex Bolt (4) until the Chair-side Adapter is secure in the Handpiece Holder.



Handpiece Holder Chair-side Adapter Assembly




Do not twist or turn the Handpiece Cable connector. The connector is keyed and mates straight on with Camera Handpiece connector.

Whether detaching just the Camera Handpiece from the Handpiece Cable and using with additional Handpiece Cables installed on other computers or bringing the complete system (handpiece and cable), the Astra System is easily shared among patient operatories.


Refer to the corresponding illustrations and perform the following when transporting only the Camera Handpiece.

1. Retract collar of quick disconnect connector and simultaneously pull straight from the Handpiece connector.
2. Release connector quick disconnect collar.
3. Store the Handpiece Cable by placing cable slide holder into the Handpiece Holder.
4. Transport Handpiece to next operatory.
5. Connect Handpiece at new location by aligning the connector key with the keyway of the mating Handpiece Cable connector.
6. Insert straight into Handpiece connector until it securely snaps into place.

When transporting the complete Astra System (handpiece and cable) among patient operatories, simply disconnect the USB connector end of the Handpiece Cable, transport the system and reconnect to the new computer. Make sure that the computer in the new operatory is using compliant software.




**Detaching Astra Handpiece**



**Connecting Astra Handpiece**

**CAUTION:** Always handle the Astra Handpiece carefully. Do not drop and make sure to protect the optical window from scratches.



**Astra Handpiece Storage**



The Astra System should be used with a Camera Cover covering the handpiece. Make sure to use a new disposable Camera Cover for each patient.







The Camera Cover must be used for only one patient and disposed of properly in accordance with local code.

**NOTE:** For optimum image clarity, the optically clear section of the camera cover must be aligned with the camera lens by facing the camera/light source section of the handpiece down toward the paper layer.

**Disposable Camera Cover Installation** - Install a new Camera Cover on the Camera Handpiece for each patient as follows:

1. Remove the Camera Handpiece from the holder.
2. Insert the Camera Handpiece tip between the second and third layers (between the outer paper and plastic layers) with the camera/light facing towards the paper layer. See A.
3. Push the Camera Handpiece completely in until the tip is fully inserted into the narrow tip of the sheath. See B.
4. Gently squeeze the illuminated tip of the Camera Handpiece between the thumb and index finger to ensure that the optically clear area of the camera cover is flat against the lens tip. See C.
5. Peel away and discard the top plastic layer and the bottom paper layer. See D.



### **Camera Cover Installation**




Verify that the video capture software is running and correctly configured.

The Astra System is operated connected directly to a computer system running various PC USB streaming video capture software applications such as VISIX. The camera provides the image/video source while the computer and associated software are used to display and save the resultant images. Perform the following procedures to operate the Astra System when correctly connected to a computer.

1. Remove the Camera Handpiece from the handpiece holder and observe that the computer display shows the high-resolution live video images taken by the camera.
2. Place the camera lens window over area of interest and view image on computer display monitor.
3. Press and release the top keypad button on the Camera Handpiece to freeze (ie. capture) the displayed image on the computer monitor screen. Press and release the top keypad button again to unfreeze the displayed image.
4. Press and release the bottom keypad button on the Camera Handpiece to save the captured image to the computer hard drive and then to return to a live image display.
5. Repeat steps 2 through 5 as necessary.
6. When imaging an X-ray film on a light box, deactivate the LED light source by depressing both Camera Handpiece buttons simultaneously for approximately 3 seconds. Perform steps 2 through 5 to freeze and save the X-ray film image as desired.  
Depress both Camera Handpiece buttons simultaneously for approximately 3 seconds to turn ON the LED light source as desired.
7. Return the Camera Handpiece to the holder when done. Verify that the LED light source and camera turns OFF.

**USB Video Diagnostic** - If live video is not displayed on the PC monitor of the Astra System when the imaging software application is running, there could be a problem with the settings in the software application or with the handpiece. Perform the following USB Video Diagnostic procedure to check that the Astra handpiece is working correctly, the Handpiece Cable is correctly connected between the Astra System and the computer and the Astra System USB driver files are installed correctly on the PC.

1. Go to the **Start** button and select **Programs**.
2. Select the **Astra** program option and click on **Camera Demonstration Program** and observe that the Camera Demonstration Program screen opens.
3. Depending on the outcome of the diagnostic, perform one of the following:
  - a. If a live image is correctly displayed, exit from the application and make sure that the imaging software application is configured correctly.
  - b. If a live image is not displayed, exit from the application and check the handpiece, Handpiece Cable and USB driver file installation.
4. Reset the driver's video settings by performing the following:
  - a. On the Astra Demonstration Program screen, click the **Properties** button located on the lower left of the live image screen.
  - b. On the Properties screen select the **VIDEO PROC AMP** tab and then click the **DEFAULT** button. The image colors should now be correct. Click **OK** to close the Properties screen and exit from application.



## INSPECTION AND CLEANING

Perform the following inspection and cleaning procedures periodically as a preventive maintenance measure to keep the Astra System in optimal condition resulting in trouble-free operation producing crisp, clear images.



Do not attempt any internal service of Astra System components. Contact your local authorized air techniques dealer for service. Failure to heed this warning may result in equipment damage and voiding the warranty.

**Inspection** - Routinely inspect each component of the Astra System for possible defects as follows:

1. Camera Handpiece -
  - a. Check overall handpiece for chips, cracks or other irregularities.
  - b. Check the lens window for debris or spots.
  - c. Check the connector socket for damage.
2. Cables and Connectors -
  - a. Check cables for damaged or deteriorated insulation kinking or twisting.
  - b. Check connectors for loose, bent or missing pins.
  - c. Check that the quick-disconnect plug end snaps close to snugly secure the Camera Handpiece connector.



Do not spray cleaning liquids or disinfectants directly on the camera handpiece. Use care not to allow liquids to run into internal circuitry.

Do not wipe the surfaces using benzine, thinner, etc. as this may degrade the finish.  
Do not autoclave the Camera Handpiece.

**CAUTION:** Unplug the Astra System components before performing cleaning.



The Camera Cover must be used for only one patient and disposed of properly in accordance with local code.

**Cleaning** - Clean each component of the Astra System as necessary by performing the following procedures.

1. Camera Handpiece - Following every patient use, remove the Camera Cover from the Camera Handpiece and dispose in accordance with local regulations. The Camera Handpiece can be wiped clean and disinfected with any non-corrosive surface disinfectant such as Birex® se.
2. Camera Window Lens - Remove debris or spots from lens window by cleaning with alcohol and cotton swabs.



---

## ***SPECIFICATIONS***

Input Power:	USB powered device using less than 500 mA. No external power supply is required. No mains connection.
Video Outputs:	High Speed USB 2.0 Video
Sensor:	¼ Inch CCD
Format:	PAL
Pixels:	720 H X 576 V
Exposure Control:	Automatic
Illumination:	Ultra-bright white LED lamps
Image:	Normal, not mirrored
Operating Temperature:	10 to 40° C (50 to 104° F)
Storage Temperature:	0 to 70° C (32 to 158° F)
Humidity:	10 to 90% non-condensing
Classifications:	Portable, Continuous Operation, Type BF Applied Part. No mains connection.
Protection for Ingress of Water:	Ordinary (per IEC 60529)
Flammable Atmosphere:	Cannot use in the vicinity of flammable anesthetic mixtures of air, oxygen or nitrous oxide.
Device Class:	USA - 21CFR/Class I Canada - MDR/Class I, EU - 93/42/EEC/Class I

---

## ***PHYSICAL CHARACTERISTICS***

<b>Dimensions</b>	<b>Length</b>	<b>Width</b>	<b>Height</b>	<b>Weight</b>
Astra Handpiece: D8060	8.6 inches (21.84 cm)	0.83 inches (2.11 cm)	1.15 inches (2.92 cm)	2.4 oz. (68.04 g)
Handpiece Cable D4020:	10 feet (3.05 m)			

# ACCESSORIES

---

The following lists the ordering number and description for accessory components available to maintain and expand the Astra System to meet your professional needs. Contact your Air Techniques Dealer for information.

---

Order Number	Description
D8000	Complete Astra System consisting of the following:
D8060	Astra Camera Handpiece Assembly
D4020	Handpiece Cable
D4170	Handpiece Holder with Chair-side Adapter
Replacement Disposable Camera Covers for Camera Handpiece	
D4125	Box of 100
D4120	Box of 500
74500	VISIX Imaging Software, 5 Licenses

## VISIX IMAGING SOFTWARE

---

Consider VISIX, our new comprehensive digital imaging software application, providing seamless integration with our digital imaging suite; the Accent CCD Intraoral Sensor, ScanX PSP Systems and Astra Intraoral Cameras.

### Key Features.

- ☐ Easy to learn.
- ☐ Simple to use.
- ☐ Acquire, mount, view and store images with just a few mouse clicks.
- ☐ Customize images to personal preferences with our simple Set-Up 1-2-3 tool.
- ☐ Viewing - Provides automatic image mounting that can be easily customized.
- ☐ Flexibility - VISIX includes a built-in bridging solution to over 35 Practice Management Software titles.

VISIX is the superb software complement for Astra. Contact your authorized dealer for further information.

This Air Techniques equipment is warranted to be free from defects in material and workmanship from the date of installation for a period of 3 years.

Any item returned to our factory through your authorized Air Techniques dealer, will be repaired or replaced at our option at no charge provided that our inspection shall indicate it to have been defective. Dealer labor, shipping and handling charges are not covered by this warranty.

This warranty does not apply to damage due to shipping, misuse, careless handling or repairs by other than authorized service personnel. Warranty void if installed or serviced by other than your authorized Air Techniques dealer service personnel. Air Techniques, Inc. is not liable for indirect or consequential damage or loss of any nature in connection with this equipment.

This warranty is in lieu of all other warranties expressed or implied. No representative or person is authorized to assume for us any liability in connection with the sale of our equipment.

---

### ON-LINE WARRANTY REGISTRATION

---

Quickly and easily register the Astra Camera on-line. Just have your product model and serial numbers available. Then go to the Air Techniques web site, **[www.airtechniques.com](http://www.airtechniques.com)**, click the **warranty registration** link and complete the registration form. This on-line registration ensures a record for the warranty period and helps us keep you informed of product updates and other valuable information.



Air Techniques, Inc.  
1295 Walt Whitman Road  
Melville, New York 11747 USA



Medical Device Safety Service  
Schiffgraben 41  
30175 Hannover, Germany

Air Techniques is a leading manufacturer of dental products including air and vacuum systems, film processors, digital imaging systems and intra-oral digital video systems. We have been manufacturing quality products for dental professionals since 1962.

Air Techniques' products are only distributed through authorized dealers.

- ☐ A/T 2000<sup>®</sup> XR
- ☐ Acadia<sup>™</sup>
- ☐ Accent<sup>™</sup>
- ☐ AirStar<sup>®</sup>
- ☐ CleanStream<sup>™</sup>
- ☐ Guardian<sup>™</sup>
- ☐ Peri-Pro<sup>®</sup>
- ☐ Astra<sup>™</sup>
- ☐ ScanX<sup>®</sup>
- ☐ Sensa<sup>™</sup>
- ☐ Spectra<sup>™</sup>
- ☐ STS<sup>™</sup>
- ☐ VacStar<sup>™</sup>
- ☐ Visix<sup>™</sup>



1295 Walt Whitman Road, Melville, NY 11747-3062

Website: [www.allproimaging.com](http://www.allproimaging.com)

East Coast: 1-800-AIR-TECH (1-800-247-8324) • West Coast: 1-800-822-2899